

TrachVox™

TrachVox

*SPEECH AND HUMIDIFICATION
FOR SPONTANEOUSLY BREATHING
TRACHEOSTOMY PATIENTS*

SMALL, COMPACT AND LIGHTWEIGHT:

- Tracheal HME with built-in speech valve
- Integrated suction port acts as relief valve during patient coughing or occlusion
- Integrated oxygen connector



*Innovative
Technology for
Tracheostomy
Patients*

For more information about **TrachVox** and our complete line of airway products visit our website at www.pulmodyne.com or call **317-246-5505**.

Pulmodyne®



TrachVox is shown actual size



Performance:

Application:

Recommended Tidal Volume Range:

Compressible Volume:

Moisture Output:

Resistance to Flow:

Weight:

Length:

Height:

Width:

Connections:

Value:

For spontaneously breathing patients using Tracheostomy Tube in either clinical or home care.

500-1000 ml

9.5 ml dead space

20.1mg H₂O/L air @V_T 500 ml

.23 +/- .15 cm H₂O @ 30 l/m

3.0 g

22mm

28mm

38mm

15mm/F TracheostomyTube; 4mm Supplemental O₂ Tubing

Information:

Connect TrachVox at open end of Tracheostomy Tube. Check that connection is properly fixed. Connect O₂ tubing to port, if required. Should be replaced every 24 hours, or as specified.

The TrachVox contains a valve with a spring that can be easily be closed airtight with a finger during speech. After releasing the finger, the valve will automatically open.

Contraindications:




Do Not Use above recommended tidal volume range. (Increased dead space can cause CO₂ retention at the low end, and reduce humidification at the high end.)

Do Not Use with patients with heavy levels of secretions from airway, or patients suffering from dehydration.

Cautions:

- **For single use only**, do not sterilize or reuse
- Before use, check for product deformity and valve functioning by depressing the top downward and observing it return to original shape
- Store in clean, dry, dark conditions
- **Discard after use**

Procedural Sequence



Caution: Federal Law (USA) restricts this device to sale by or on the order of a licensed physician.

Single Patient Use

Rx ONLY

Latex Free

Pulmodyne[®]

... bringing change to life[®]

Pulmodyne
2055 Executive Drive
Indianapolis, IN 46241 USA

T: 317.246.5505

F: 317.246.5501

www.Pulmodyne.com

Pulmodyne and bringing change to life are registered trademarks of Pulmodyne, Inc. TrachVox is a trademark of Pulmodyne, Inc.